0	9/97 (9.3%)	35/89 (39.3%)	44/78 (56.4%)	34/57 (59.6%)	47/51 (92.2%)	44/47 (93.6%)
2	17/34 (50.0%)	23/34 (67.6%)	25/32 (78.1%)	58/62 (93.5%)	59/60 (98.3%)	43/43 (100%)
4	16/21 (76.2%)	15/20 (75.0%)	13/15 (86.7%)	55/58 (94.8%)	53/54 (98.1%)	39/40 (97.5%)
6	20/28 (71.4%)	21/27 (77.8%)	16/18 (88.9%)	27/34 (79.4%)	26/33 (78.8%)	22/27 (81.5%)
8	13/23 (56.5%)	13/22 (59.1%)	4/11 (36.4%)	12/17 (70.6%)	12/17 (70.6%)	7/11 (63.6%)
10	5/9 (55.6%)	5/8 (62.5%)		4/9 (44.4%)	3/7 (42.9%)	

1 Test= one (1) test performed on the noted days after first PCR positive test result. Day 0 is the first day of documented infection with SARS-CoV-2.
2 Tests= two (2) tests performed an average of 48 hours apart. The first test performed on the indicated day and the second test performed 48 hours later.
3 Tests= three (3) tests performance an average of 48 hours apart. The first test performed on the indicated day, the second test performed 48 hours later, and a final test performed 48 hours after the second test.

Cross Reactivity:

Cross reactivity with the following organisms has been studied. Samples positive for the following organisms were found negative when tested with the Rapid Response™ COVID-19 Antigen Rapid Test Device.

HCoV-HKU1	Influenza A (H5N1)	Coxsackie virus A16 (Enterovirus)
HCoV-OC43	Influenza A (H7N9)	Norovirus
HCoV-NL63	Influenza A (H7N7)	Mump virus
HCoV-229E	Influenza B Victoria lineage	<i>Legionella pneumophila</i>
Measles virus	Influenza B Yamagata lineage	<i>Mycoplasma pneumoniae</i>
<i>Streptococcus pneumoniae</i>	Respiratory syncytial virus	<i>Chlamydia pneumoniae</i>
Epstein-Barr virus	Adenovirus	<i>Streptococcus pyogenes</i>
Bordetella Para pertussis	Parainfluenza 1/2/3 virus	<i>Streptococcus agalactiae</i>
Influenza A (H1N1) pdm09	Human metapneumovirus	Group C <i>Streptococcus</i>
Influenza A (H3N2)	Rhinovirus	<i>Staphylococcus aureus</i>
MERS-coronavirus	Parainfluenza virus 4	<i>Staphylococcus epidermidis</i>
<i>Haemophilus influenzae</i>	<i>Candida albicans</i>	

Microbial Interference Study:

Potential microbial interference was evaluated to demonstrate that false negatives will not occur when SARS-CoV-2 is present in a specimen with other microorganisms. Low concentration of SARS-CoV-2 (3 X LOD) was spiked into the higher concentrations of interfering organism, and it was found that there is no microbial interference for following organisms.

HCoV-HKU1	Influenza A (H5N1)	Coxsackie virus A16 (Enterovirus)
HCoV-OC43	Influenza A (H7N9)	<i>Haemophilus influenzae</i>
HCoV-NL63	Influenza A (H7N7)	<i>Candida albicans</i>
HCoV-229E	Influenza B Victoria lineage	<i>Mycobacterium tuberculosis</i>
Measles virus	Influenza B Yamagata lineage	Norovirus
<i>Streptococcus pneumoniae</i>	Respiratory syncytial virus	Mump virus
Epstein-Barr virus	Adenovirus	<i>Legionella pneumophila</i>
Bordetella Para pertussis	Parainfluenza 1/2/3 virus	<i>Mycoplasma pneumoniae</i>
Influenza A (H1N1) pdm09	Human metapneumovirus	<i>Chlamydia pneumoniae</i>
Influenza A (H3N2)	Rhinovirus	<i>Streptococcus pyogenes</i>
Group C <i>Streptococcus</i>	<i>Staphylococcus aureus</i>	<i>Streptococcus agalactiae</i>
SARS-coronavirus	Parainfluenza virus 4	<i>Staphylococcus epidermidis</i>
MERS-coronavirus	<i>Candida albicans</i>	<i>Haemophilus influenzae</i>
Pooled human nasal wash – representative of normal respiratory microbial flora		

Interfering Substances


The following substances, naturally present in respiratory specimens or that may be artificially introduced into the respiratory tract, were evaluated at the concentrations listed below. None of them were found to affect test performance of the Rapid Response™ COVID-19 Antigen Rapid Test Device.

Substance	Concentration	Substance	Concentration
3 OTC nasal sprays	10%	Guaiacol glyceryl ether	20 mg/ml
3 OTC mouthwashes	10%	Mucin	1%
3 OTC throat drops	10%	Mupirocin	250 µg/ml
4-acetamidophenol	10 mg/ml	Oxymetazoline	10 mg/ml
Acetylsalicylic acid	20 mg/ml	Phenylephrine	10 mg/ml
Albuterol	20 mg/ml	Phenylpropanolamine	20 mg/ml
Chlorpheniramine	5 mg/ml	Relenza® (zanamivir)	20 mg/ml
Dexamethasone	5 mg/ml	Rimantadine	500 ng/ml
Dextromethorphan	10 mg/ml	Tamiflu® (oseltamivir)	100 mg/ml
Diphenhydramine	5 mg/ml	Tobramycin	40 mg/ml
Doxylamine succinate	1 mg/ml	Triamcinolone	14 mg/ml
Flunisolide	3 mg/ml		

High Dose Hook Effect

No high dose hook effect was observed when tested with up to a concentration of 1 x 10^{6.4} TCID₅₀/mL of heat inactivated SARS-CoV-2 virus with the Rapid Response™ COVID-19 Rapid Test Device.

GLOSSARY OF SYMBOLS

	Consult instructions for use		Tests per Kit		Catalogue number
	Store between 2°C to 30°C		Use by date		Do Not Reuse
	In vitro diagnostic medical device		Lot Number		Authorized Representative

BTNX, Inc.
722 Rosebank Rd
Pickering, ON, L1W 4B2, Canada
Technical Support: 1-888-339-9964

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**MDSS GmbH**